

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Denhardt´s Solution, 50x Konzentrat	Kat.Nr./Cat.No. : 39603
Denhardt´s solution, 50x concentrate molecular biology grade	Lot/Contr.No. : P120473

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	klare farblose bis leicht gelbliche Lösung clear colorless to slightly yellowish solution	entspricht corresponds
pH		7.5
Zusammensetzung composition		
Albumin Fraktion V	0.1g/ 10 ml	entspricht corresponds
Polyvinylpyrrolidone	0.1g/ 10 ml	entspricht corresponds
Ficoll [®] 400	0.1g/ 10 ml	entspricht corresponds
DNase/ RNase Aktivität DNase/ RNase activity	nicht nachweisbar none detected	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		11/ 2014
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 22.11.2012

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.