

**A N A L Y S E N Z E R T I F I K A T**  
CERTIFICATE OF ANALYSIS

<b>Denhardt´s Solution, 50x Konzentrat</b>	<b>Kat.Nr./Cat.No. : 39603</b>
Denhardt´s solution, 50x concentrate	<b>Lot/Contr.No. : P220009</b>
molecular biology grade	

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Aussehen</b> appearance	klare farblose bis leicht gelbliche Lösung clear colorless to slightly yellowish solution	entspricht corresponds
<b>pH</b>		7.0
<b>Zusammensetzung</b> composition		
Albumin Fraktion V	0.1g/ 10 ml	entspricht corresponds
Polyvinylpyrrolidone	0.1g/ 10 ml	entspricht corresponds
Ficoll <sup>®</sup> 400	0.1g/ 10 ml	entspricht corresponds
<b>DNase/ RNase Aktivität</b> DNase/ RNase activity	nicht nachweisbar none detected	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life		13.01.2025
<b>Lagerung (°C)</b> storage		-15 bis/ to -25

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 20.11.2023**

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.