

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

EDTA-Na₂-Lösung, 0.5M	Kat.Nr./Cat.No. : 39761
Ethylenediamine tetraacetic acid Na ₂ -salt solution, 0.5M molecular biology grade	Lot/Contr.No. : P220150

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	klare, farblose Lösung clear, colorless solution	entspricht corresponds
pH	7.9 – 8.1	8.1
Zusammensetzung composition		
EDTA · Na ₂ · 2H ₂ O (#39760)	186.1 g/l	entspricht corresponds
NaOH (10M)	ca. 50 ml/l	entspricht corresponds
DNase/ RNase Aktivität DNase/ RNase activity	nicht nachweisbar none detected	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		14.04.2024
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 27.04.2022

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.