

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Tris(hydroxymethyl)aminomethan-hydrochlorid	Kat.Nr./Cat.No. :	39787
TRIS-HCL molecular biology grade	Lot/Contr.No. :	160916

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₄ H ₁₁ NO ₃ ·HCl	
Molekulargewicht molecular weight		157.60	
Aussehen appearance		weißes kristallines Pulver white crystalline powder	entspricht corresponds
Gehalt (%) assay	Titration	min. 99.0	100.5
Löslichkeit solubility	10 % in H ₂ O	klare, farblose Lösung clear and colorless solution	entspricht corresponds
pH	10 % in H ₂ O	3.5 - 5.0	4.2
A 1cm/ 230 nm	10 % in H ₂ O	max. 0.1	0.046
A 1cm/ 260 nm	10 % in H ₂ O	max. 0.05	0.028
A 1cm/ 280 nm	10 % in H ₂ O	max. 0.03	0.027
Schmelzpunkt (°C) melting point		148 – 153	152
Schwermetalle (ppm) heavy metals	als / as Pb	max. 10	< 5
Identität identity	IR	entspricht corresponds	entspricht corresponds
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			22.09.2021
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 27.09.2016

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.