

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Resazurin Zell-Viabilitäts-Assay	Kat.Nr./Cat.No. :	39905
Resazurin Cell Viability Assay	Lot/Contr.No. :	180595

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	dunkelblaue Lösung dark-blue solution	entspricht corresponds
Anwendung application	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		30.04.2020
Lagerung (°C) storage		-15 bis/ to -25

SERVA Electrophoresis GmbH

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 29.05.2018

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.