

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Resazurin Zell-Viabilitäts-Assay	Kat.Nr./Cat.No. :	39905
Resazurin Cell Viability Assay	Lot/Contr.No. :	220435

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	dunkelblaue Lösung dark-blue solution	entspricht corresponds
Anwendung application	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		30.04.2024
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control**Druckdatum/ printing date:** 03.05.2022

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.