

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Resazurin Zell-Viabilitäts-Assay Resazurin Cell Viability Assay	Kat.Nr./Cat.No. : 39905 Lot/Contr.No. : 2209875
---	--

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	dunkelblaue Lösung dark-blue solution	entspricht corresponds
Anwendung application	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		30.09.2024
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 13.10.2022

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.