

**A N A L Y S E N Z E R T I F I K A T**  
**C E R T I F I C A T E O F A N A L Y S I S**

<b>DOWEX® 1X2</b>	<b>Kat.Nr./Cat.No. : 41010</b>
(50 – 100 mesh) pract.	<b>Lot/Contr.No. : 150944</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		reinweißes bis elfenbeinfarbenes Korn off-white to ivory- colored particles	entspricht corresponds
<b>Trocknungsverlust (%)</b> loss on drying	3 h, 105 °C	65.0 – 75.0	67.1
<b>Gegenion</b> counter-ion		Cl <sup>-</sup>	entspricht corresponds
<b>Kapazität (eq/l)</b> capacity		min. 0.7	0.7
<b>Mindesthaltbarkeit</b> minimum shelf life			06.08.2019
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Dieser Artikel kann produktionsbedingt andersfarbige Partikel enthalten.**  
**Diese haben keine Auswirkungen auf die Eigenschaften und die Funktionalität des Produktes.**  
 Due to the production process this article can contain differently coloured particles.  
 This does not affect the quality characteristics and the functionality of the product.

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
 We do not guarantee that the product can be used for a special application.  
 This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 05.09.2017**

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.