

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>DOWEX<sup>®</sup> 1X2</b> , p.A. ( 50 – 100 mesh) analytical grade	<b>Kat.Nr./Cat.No. :</b> 41011 <b>Lot/Contr.No. :</b> P080166
--	--

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		reinweißes bis elfenbeinfarbenes Korn off-white to ivory-colored particles	entspricht corresponds
<b>Kapazität (eq/l)</b> capacity		min. 0.7	0.8
<b>Trocknungsverlust (%)</b> loss on drying	3h, 105°C	65.0 - 75.0	70.6
<b>Gegenion</b> counter-ion		Cl <sup>-</sup>	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			05/ 2016
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 02.06.2015

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.