

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Dowex® 1X2 p.A. (200 – 400 mesh) analytical grade	Kat.Nr./Cat.No. : 41031 Lot/Contr.No. : P150229
---	--

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weißes bis hellelfenbein- farbenes Korn white to light-ivory colored particles	entspricht corresponds
Trocknungsverlust (%) loss on drying	105 °C, 3 h	70.0 – 80.0	78.7
Gegenion counter ion		Cl ⁻	entspricht corresponds
Kapazität (eq/l) capacity		min. 0.6	0.6
Mindesthaltbarkeit minimum shelf life			06/ 2020
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 30.06.2015

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.