

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

DOWEX® 1X8. (20 – 50 mesh), pract.	Kat.Nr./Cat.No. : 41080
	Lot/Contr.No. : 180580

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		hellelfenbeinfarbenes bis beiges Korn light ivory-coloured to beige bead	entspricht corresponds
Trocknungsverlust (%) loss on drying	105 °C, 3 h	45.0 – 60.0	52.1
Gegenion counter-ion		Cl ⁻	entspricht corresponds
Kapazität (eq/l) capacity		min. 1.3	1.38
Mindesthaltbarkeit minimum shelf life			24.05.2023
Lagerung (°C) storage			+15 bis/ to +30

Dieser Artikel kann produktionsbedingt andersfarbige Partikel enthalten.
Diese haben keine Auswirkungen auf die Eigenschaften und die Funktionalität des Produktes.
Due to the production process this article can contain differently coloured particles.
This does not affect the quality characteristics and the functionality of the product.

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 30.05.2018

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.