

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>DOWEX® 1X8</b> , p.A. (20 – 50 mesh) analytical grade	<b>Kat.Nr./Cat.No. :</b> 41081 <b>Lot/Contr.No. :</b> P140190
-------------------------------------------------------------	------------------------------------------------------------------

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		sandgelbes, transparentes Korn sandyellow-colored, transparent particles	entspricht corresponds
<b>Trocknungsverlust (%)</b> loss on drying	3 h, 105 °C	50.0 – 60.0	51.5
<b>Gegenion</b> counter-ion		Cl <sup>-</sup>	entspricht corresponds
<b>Kapazität (eq/l)</b> capacity		min. 1.3	1.3
<b>Mindesthaltbarkeit</b> minimum shelf life			04/ 2019
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Das Produkt kann schwarze Partikel enthalten. Dabei handelt es sich um inaktive Copolymere, die die spezifizierten Eigenschaften wie Kapazität, Wassergehalt oder die Abgabe von organischen Substanzen vom Harz nicht beeinträchtigen.**

The product can contain black particles. These are inactive copolymers which do not impair the specified characteristics as capacity, water content or the release of organic substances of the resin.

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**

**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 19.05.2014

Dipl.-Ing. (FH) Bernhard Göckel

Jakob Arndt

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**

This report has been computer generated and does not contain a signature.