

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

| | |
|---|--|
| DOWEX® 1X8 , p.A. (20 – 50 mesh) analytical grade | Kat.Nr./Cat.No. : 41081 Lot/Contr.No. : P140351 |
|---|--|

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|---|-------------------|---|---------------------------|
| Aussehen appearance | | sandgelbes, transparentes Korn sandyellow-colored, transparent particles | entspricht corresponds |
| Trocknungsverlust (%) loss on drying | 3 h, 105 °C | 50.0 – 60.0 | 54.1 |
| Gegenion counter-ion | | Cl ⁻ | entspricht corresponds |
| Kapazität (eq/l) capacity | | min. 1.3 | 1.3 |
| Mindesthaltbarkeit minimum shelf life | | | 07/ 2019 |
| Lagerung (°C) storage | | | +15 bis/ to +30 |

Das Produkt kann schwarze Partikel enthalten. Dabei handelt es sich um inaktive Copolymere, die die spezifizierten Eigenschaften wie Kapazität, Wassergehalt oder die Abgabe von organischen Substanzen vom Harz nicht beeinträchtigen.

The product can contain black particles. These are inactive copolymeres which do not impair the specified characteristics as capacity, water content or the release of organic substances of the resin.

| |
|---|
| <p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p> |
|---|

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 30.07.2014

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.