

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

DOWEX â 1X8 (100 – 200 mesh) , pract.	Kat.Nr./Cat.No. : 41100
	Lot/Contr.No. : 170984

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		elfenbeinfarbenes Korn ivory colored particles	entspricht corresponds
Gegenion counter ion		Cl ⁻	entspricht corresponds
Trocknungsverlust (%) loss on drying	3 h; 105 °C	39.0 – 45.0	42.3
Kapazität (eq/ l) capacity		min. 1.2	1.5
Mindesthaltbarkeit minimum shelf life			03.11.2022
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 08.11.2017

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.