

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

DOWEX[®] 50 WX2 p.A. (100 – 200 mesh), analytical grade	Kat.Nr./Cat.No. : 41521 Lot/Contr.No. : P160340
--	--

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		elfenbeinfarbenes Korn ivory-colored particles	entspricht corresponds
Trocknungsverlust (%) loss on drying	3 h, 105 °C	74.0 – 82.0	80.0
Gegenion counter ion		Na ⁺	entspricht corresponds
Kapazität (eq/l) capacity		min. 0.6	0.7
Mindesthaltbarkeit minimum shelf life			19.08.2021
Lagerung (°C) storage			+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 23.08.2016

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.