

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

DOWEX® 50 WX2 p.A. (100 – 200 mesh), analytical grade	Kat.Nr./Cat.No. : 41521 Lot/Contr.No. : P200029
---	--

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		elfenbeinfarbenes Korn ivory-colored particles	entspricht corresponds
Trocknungsverlust (%) loss on drying	3 h, 105°C	74.0 – 82.0	81.6
Gegenion counter ion		Na ⁺	entspricht corresponds
Kapazität (eq/l) capacity		min. 0.6	0.7
Mindesthaltbarkeit minimum shelf life			23.01.2025
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 24.01.2020

Christian Monsler

Marina Konter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.