

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

| | |
|---|---------------------------------------|
| DOWEX® 50 WX8 (100 – 200 mesh) pract. | Kat.Nr./Cat.No. : 41620 |
| | Lot/Contr.No. : 161177 |

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|---|-------------------|--|---------------------------|
| Aussehen appearance | | beige-grünes Korn beige-green-colored particles | entspricht corresponds |
| Trocknungsverlust (%) loss on drying | 3 h, 105 °C | 50.0 – 58.0 | 52.1 |
| Gegenion counter-ion | | H ⁺ | entspricht corresponds |
| Kapazität (eq/l) capacity | | min. 1.7 | 2.0 |
| Mindesthaltbarkeit minimum shelf life | | | 27.07.2019 |
| Lagerung (°C) storage | | | +15 bis/ to +30 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 16.01.2017

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.