

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

DOWEX® 50 WX8 p.A.	Kat.Nr./Cat.No. : 41621
(100 – 200 mesh) analytical grade	Lot/Contr.No. : P140149

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		beiges Korn beige-colored particles	entspricht corresponds
Trocknungsverlust (%) loss on drying	3 h, 105 °C	45.0 – 55.0	46.2
Gegenion counter-ion		Na ⁺	entspricht corresponds
Kapazität (eq/l) capacity		min. 1.7	2.1
Mindesthaltbarkeit minimum shelf life			03/ 2019
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 21.03.2014

Dipl.-Ing. (FH) Bernhard Göckel

Jakob Arndt

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.