

**A N A L Y S E N Z E R T I F I K A T**  
CERTIFICATE OF ANALYSIS

<b>SERVA Tris-Glycin nativer Elektrophorese Puffer (10x)</b>	<b>Kat.Nr./ Cat.No. : 42530</b>
<b>SERVA Tris-Glycine native electrophoresis buffer (10x)</b>	<b>Lot/ Contr.No. : P200268</b>

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Aussehen</b> appearance	klare farblose Lösung clear colorless solution	entspricht corresponds
<b>pH</b>	8.4 – 8.8	8.7
<b>Zusammensetzung</b> composition		
TRIS (# 37180)	0.25 M	entspricht corresponds
Glycin/ Glycine (# 23390)	1.92 M	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life		31.01.2022
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 31.07.2020**

Christian Monsler

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.