

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

SERVA Blue G Lösung für BN, 1%	Kat.Nr./Cat.No. :	42538
SERVA Blue G solution for BN, 1%	Lot/Contr.No. :	P160015

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	dunkelblaue Lösung darkblue solution	entspricht corresponds
Zusammensetzung composition		
Serva Blue G250 (#35050)	1 %	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		14.01.2020
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 06.12.2017

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.