

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

SERVA Blue G Lösung für BN, 1%	Kat.Nr./Cat.No. : 42538
SERVA Blue G solution for BN, 1%	Lot/Contr.No. : P190013

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	dunkelblaue Lösung darkblue solution	entspricht corresponds
Zusammensetzung composition		
Serva Blue G250 (#35050)	1 %	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		10.01.2023
Lagerung (°C) storage		+15 bis/ to +30

SERVA Electrophoresis GmbH

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 10.11.2020

Christian Monsler

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.