

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS

SERVAGel(TM) IEF Running Buffer Kit

Kat.Nr./Cat.No. : 42539

Lot/Contr.No. : P250314

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	Für die isoelektrische Fokussierung von Proteinen For the isoelectric focusing of proteins	entspricht corresponds
Inhalt content	SERVA IEF Anode Buffer, powder SERVA IEF Cathode Buffer, powder	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		23.05.2027
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 25.09.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.