

ANALYSENZERTIFIKAT
 CERTIFICATE OF ANALYSIS

SERVA Tris-SDS Probenpuffer pH 8.6 (10X)	Kat.Nr./Cat.No. : 42546
SERVA Tris-SDS sample buffer pH 8.6 (10X)	Lot/Contr.No. : P190113

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		blau-rote Lösung blue-red solution	entspricht corresponds
pH	1x konzentriert 1x concentrated	8.5 – 8.8	8.5
Zusammensetzung composition		entspricht corresponds	entspricht corresponds
Glycerin glycerol		25%	
TRIS		290 mmol/ l	
TRIS-HCl		210 mmol/ l	
Natriumdodecylsulfat sodium dodecyl sulfate		10%	
Bromphenolblau bromophenol blue		0.1%	
Mindesthaltbarkeit minimum shelf life			13.09.2024
Lagerung (°C) storage			+2 bis/ to +8

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control
Druckdatum/ printing date: 20.07.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.