

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

---

|                         |                          |                |
|-------------------------|--------------------------|----------------|
| <b>TE Puffer (100x)</b> | <b>Kat.Nr./Cat.No. :</b> | <b>42554</b>   |
| TE buffer (100x)        | <b>Lot/Contr.No. :</b>   | <b>P220034</b> |

---

| <b>Parameter</b><br>parameter                     | <b>Spezifikation</b><br>specification             | <b>Ergebnis</b><br>result |
|---|---|---------------------------|
| <b>Aussehen</b><br>appearance                     | klare farblose Lösung<br>clear colorless solution | entspricht<br>corresponds |
| <b>pH</b>   | 7.8 – 8.2   | 8.2                       |
| <b>Zusammensetzung</b><br>composition             |   |                           |
| TRIS (# 37180)                                    | 1.0 M   | entspricht<br>corresponds |
| EDTA-Na <sub>2</sub> ·2H <sub>2</sub> O (# 11280) | 0.1 M   | entspricht<br>corresponds |
| <b>Mindesthaltbarkeit</b><br>minimum shelf life   |   | 04.02.2024                |
| <b>Lagerung (°C)</b><br>storage                   |   | +15 bis/ to +30           |

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

---

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control****Druckdatum/ printing date:** 09.02.2022

Christian Monsler

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.