

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

|                                          |                                |
|------------------------------------------|--------------------------------|
| <b>SSC-Puffer (20x)</b>                  | <b>Kat.Nr./Cat.No. : 42555</b> |
| SSC-buffer (20x) molecular biology grade | <b>Lot/Contr.No. : P210118</b> |

| <b>Parameter</b><br>parameter                                                           | <b>Spezifikation</b><br>specification             | <b>Ergebnis</b><br>result |
|-----------------------------------------------------------------------------------------|---------------------------------------------------|---------------------------|
| <b>Aussehen</b><br>appearance                                                           | klare farblose Lösung<br>clear colorless solution | entspricht<br>corresponds |
| <b>pH</b>                                                                               | 6.8 - 7.2                                         | 6.9                       |
| <b>Zusammensetzung</b><br>composition                                                   |                                                   |                           |
| NaCl (# 30183)                                                                          | 3.0 M                                             | entspricht<br>corresponds |
| tri-Natriumcitrat·2H <sub>2</sub> O<br>tri-Sodiumcitrate·2H <sub>2</sub> O<br>(# 38642) | 0.3 M                                             | entspricht<br>corresponds |
| <b>DNase/ RNase Aktivität</b><br>DNase/ RNase activity                                  | nicht nachweisbar<br>none detected                | entspricht<br>corresponds |
| <b>Mindesthaltbarkeit</b><br>minimum shelf life                                         |                                                   | 14.10.2024                |
| <b>Lagerung (°C)</b><br>storage                                                         |                                                   | +15 bis/ to +30           |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.  
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 21.08.2023**

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.  
This report has been computer generated and does not contain a signature.