

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>TBE Puffer 10x</b>	<b>Kat.Nr./Cat.No. : 42557</b>
TBE Buffer 10x	<b>Lot/Contr.No. : P120387</b>

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Aussehen</b> appearance	klare farblose Lösung clear colorless solution	entspricht corresponds
<b>pH</b>	8.1 – 8.5	8.4
<b>Zusammensetzung</b> composition		
TRIS (# 37180)	0.89 M	entspricht corresponds
EDTA-Na <sub>2</sub> ·2H <sub>2</sub> O (# 11280)	0.02 M	entspricht corresponds
Borsäure / Boric acid (# 15165)	0.89 M	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life		09/ 2014
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 27.09.2012**

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Steinmann

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.