

CERTIFICATE OF ANALYSIS

Semi-dry Blotting Buffer Kit	Cat.No. :	42559
for Western Blotting	Contr.No. :	P130047

Kit Components

Buffer I	conc. anode buffer	0.3 M Tris and 20 % methanol in aqueous solution
Buffer II	diluted anode buffer	0.03 M Tris and 20 % methanol in aqueous solution
Buffer III	cathode buffer	0.025 M Tris/HCl pH 9.4, 0.04 M 6-aminocaproic acid and 20 % methanol in aqueous solution

Application

Western Blotting

Minimum shelf life

02/ 2018

Storage (°C)

+15 to +30

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 19.02.2013

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.