

CERTIFICATE OF ANALYSIS

Semi-dry Blotting Buffer Kit for Western Blotting	Cat.No. :	42559
	Contr.No. :	P140440

Kit Components

Buffer I	conc. anode buffer	0.3 M Tris and 20 % methanol in aqueous solution
Buffer II	diluted anode buffer	0.03 M Tris and 20 % methanol in aqueous solution
Buffer III	cathode buffer	0.025 M Tris/HCl pH 9.4, 0.04 M 6-aminocaproic acid and 20 % methanol in aqueous solution

Application

Western Blotting

Minimum shelf life

09/ 2019

Storage (°C)

+15 to +30

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 15.10.2014

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

This report has been computer generated and does not contain a signature.