

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>SERVA Tris-Tricin/SDS Elektrophorese Puffer (20x)</b>	<b>Kat.Nr./Cat.No. :</b>	<b>42560</b>
<b>SERVA Tris-Tricine/SDS electrophoresis buffer (20x)</b>	<b>Lot/Contr.No. :</b>	<b>P110302</b>

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Aussehen</b> appearance	klare farblose Lösung clear colorless solution	entspricht corresponds
<b>pH</b>	8.3 – 8.7	8.6
<b>Zusammensetzung</b> composition		
TRIS (# 37180)	1.2 M	entspricht corresponds
Tricin/ Tricine (# 37195)	800 mM	entspricht corresponds
SDS (# 20763)	2 %	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life		07/ 2016
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 29.05.2015**

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.