

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS

BlueBlock PF (10x)	Kat.Nr./Cat.No. : 42591
BlueBlock PF (10x)	Lot/Contr.No. : P250247

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		leicht gelbliche Lösung slightly yellowish solution	entspricht corresponds
pH	1:10 verdünnt/ diluted in H ₂ O	7.1 – 7.5	7.2
Anwendung application		entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			31.07.2027
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 22.08.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.