

**ANALYSENZERTIFIKAT**  
 CERTIFICATE OF ANALYSIS

<b>PBS Puffer (10x), steril</b>	<b>Kat.Nr./Cat.No. : 42595</b>
PBS Buffer (10x), sterile	<b>Lot/Contr.No. : P190173</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		klare, farblose Lösung clear, colourless solution	entspricht corresponds
<b>pH</b>	1:10 verdünnt 1:10 diluted	7.2 – 7.6	7.4
<b>Steriltest</b> sterile test		kein Wachstum no growth	entspricht corresponds

**Zusammensetzung**  
 composition

NaCl (# 30183)	1.37 M
KCl (# 26868)	27 mM
Na <sub>2</sub> HPO <sub>4</sub> •2H <sub>2</sub> O (# 30200)	100 mM
KH <sub>2</sub> PO <sub>4</sub> (# 26870)	20 mM

<b>Mindesthaltbarkeit</b> minimum shelf life	02.05.2021
---	------------

<b>Lagerung (°C)</b> storage	+15 bis/ to +30
---------------------------------	-----------------

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b>  <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b>          We do not guarantee that the product can be used for a special application.          This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**
**Druckdatum/ printing date: 13.05.2019**

Christian Monsler

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.