

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

BlueClear SB, für das Strippen von Antikörpern	Kat.Nr./Cat.No. : 42599
BlueClear SB, for antibody stripping	Lot/Contr.No. : P220215

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, farblose Lösung clear, colourless solution	entspricht corresponds
pH		2.6 – 2.8	2.7
Mindesthaltbarkeit minimum shelf life			21.06.2023
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control**Druckdatum/ printing date: 04.07.2022**

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.