

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Xpress NC Blotting Kit for Western Blotting	Kat.Nr./Cat.No.: 42663 Lot/Contr.No.: P150424
---	--

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
------------------------	-------------------	--------------------------------	--------------------

Beschreibung description	Kit für schnelles Semi-Dry Western Blotting von 10 vertikalen mini SDS PAGE Gelen. Kit for fast Semi-Dry Western Blotting of 10 vertical mini SDS PAGE gels.		
------------------------------------	---	--	--

Kitkomponenten
kit-components

250 ml SERVA Xpress Blotting Buffer (10x)

10 Connection Paper (80 x 85 mm)

20 Blotting Fleece sheets (80 x 85 mm)

10 NC 2 Supported Nitrocellulose Membrane
(Pore size 0.2 µm, 80 x 85 mm)

Mindesthaltbarkeit minimum shelf life	27.08.2017
---	------------

Lagerung (°C) storage	+15 bis/ to +30
---------------------------------	-----------------

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 08.09.2016

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.