

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Xpress PVDF Blotting Kit	Kat.Nr./Cat.No. : 42664
for Western Blotting	Lot/Contr.No. : P180216

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
------------------------	-------------------	--------------------------------	--------------------

Beschreibung description	Kit für schnelles Semi-Dry Western Blotting von 10 vertikalen mini SDS PAGE Gelen. Kit for fast Semi-Dry Western Blotting of 10 vertical mini SDS PAGE gels.
------------------------------------	---

Kitkomponenten kit-components	SERVA Xpress Blotting Buffer (10x) Connection Paper (80 x 85 mm) Blotting Fleece sheets (80 x 85 mm) PVDF 0.2 transfer membrane (Pore size: 0.2 µm, 80 x 85 mm)
---	---

Mindesthaltbarkeit minimum shelf life	05.01.2021
---	------------

Lagerung (°C) storage	+15 bis/ to +30
---------------------------------	-----------------

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 03.01.2020

Christian Monsler

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.