

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Xpress PVDF Blotting Kit for Western Blotting	Kat.Nr./Cat.No. : 42664 Lot/Contr.No. : P200250
---	--

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
-------------------------------	--------------------------	---------------------------------------	---------------------------

Beschreibung description	Kit für schnelles Semi-Dry Western Blotting von 10 vertikalen mini SDS PAGE Gelen. Kit for fast Semi-Dry Western Blotting of 10 vertical mini SDS PAGE gels.
------------------------------------	---

Kitkomponenten kit-components	<p>SERVA Xpress Blotting Buffer (10x)</p> <p>Connection Paper (80 x 85 mm)</p> <p>Blotting Fleece sheets (80 x 85 mm)</p> <p>PVDF 0.2 transfer membrane (Pore size: 0.2 µm, 80 x 85 mm)</p>
---	---

Mindesthaltbarkeit minimum shelf life	21.03.2022
---	------------

Lagerung (°C) storage	+15 bis/ to +30
---------------------------------	-----------------

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 31.07.2020

Dr. Judith Koch

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.