

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Xpress PVDF Blotting Kit
for Western Blotting

Kat.Nr./Cat.No. : 42664
Lot/Contr.No. : P240096

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	Kit für schnelles Semi-Dry Western Blotting von 10 vertikalen mini SDS PAGE Gelen. Kit for fast Semi-Dry Western Blotting of 10 vertical mini SDS PAGE gels.		
Kitkomponenten kit-components	SERVA Xpress Blotting Buffer (10x) Connection Paper (80 x 85 mm) Blotting Fleece sheets (80 x 85 mm) PVDF 0.2 transfer membrane (Pore size: 0.2 µm, 80 x 85 mm)		
Mindesthaltbarkeit minimum shelf life	28.03.2025		
Lagerung (°C) storage	+15 bis/ to +30		

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 18.03.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.