

**A N A L Y S E N Z E R T I F I K A T**  
 C E R T I F I C A T E O F A N A L Y S I S

<b>Servalyt® 5 - 6</b>	<b>Kat.Nr./Cat.No. :</b>	<b>42924</b>
	<b>Lot/Contr.No. :</b>	<b>P190260</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Beschreibung</b> description	Klare Lösung von Polyamino-polysulfon-polycarboxylsäuren, 40% (w/v) in Wasser, sterilfiltriert Clear solution of polyamino-polysulfonic-polycarboxylic acids, 40% (w/v) in water, sterile filtered		
<b>pH</b>	2% in H <sub>2</sub> O, RT	5.53 ± 0.5	5.6
<b>Leitfähigkeit ( mS/cm )</b> conductivity	2% in H <sub>2</sub> O, RT	0.100 – 0.290	0.176
<b>pH-Bereich</b> pH-range	IEF <sup>(1)</sup> , 150 µm-Gel, 5°C, 5 000 Vh	4.8 ± 0.3 - 6.3 ± 0.3	4.7 - 6.0
<b>ausgebildeter pH-Gradient (%)</b> pH-gradient formed	IEF, 150 µm-Gel, 5°C, 5 000 Vh	min. 60	69
<b>Mindesthaltbarkeit</b> minimum shelf life			02.02.2024
<b>Lagerung (°C)</b> storage			+2 bis/ to +8

<sup>(1)</sup> Gelelektrophorese wird auf Standard-PAG-Gelen durchgeführt.  
 Gel electrophoresis is performed using a standard-PAG-gel.

( 5% T, 3% C, 4% Servalyt<sup>®</sup>, 0.1% (w/v) APS, 0.078% (w/v) TEMED,  
 standard electrode solution # 42984 and # 42986)

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
 We do not guarantee that the product can be used for a special application.  
 This document does not release you from performing the standard control upon receipt of incoming goods.

<b>SERVA Electrophoresis GmbH</b>	<b>Druckdatum/ printing date: 16.12.2020</b>
<b>Qualitätskontrolle/ Quality Control</b>	

Christian Monsler

Marina Konter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.