

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS

SERVALYT® 6 - 7
Kat.Nr./Cat.No. : 42925
Lot/Contr.No. : P090319

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	Klare Lösung von Polyamino-polysulfon-polycarboxylsäuren, 40% (w/v) in Wasser, sterilfiltriert Clear solution of polyamino-polysulfonic-polycarboxylic acids, 40% (w/v) in water, sterile filtered		
pH	2% in H ₂ O, RT	6.7 ± 0.5	6.6
Leitfähigkeit (mS/cm) conductivity	2% in H ₂ O, RT	0.100 – 0.450	0.370
pH-Bereich pH-range	IEF ⁽¹⁾ , 150 µm-Gel, 5°C, 5 000 Vh	5.6 ± 0.3 - 7.4 ± 0.3	5.9 – 7.6
ausgebildeter pH-Gradient (%) pH-gradient formed	IEF, 150 µm-Gel, 5°C, 5 000 Vh	min. 70	70.6
Mindesthaltbarkeit minimum shelf life			06.10.2022
Lagerung (°C) storage			+2 bis/ to +8
Hinweis note	aus produktionstechnischen Gründen ist die Farbe dieses Servalyt® 6-7 um einige Nuancen dunkler. Die Qualität wird dadurch nicht beeinflusst. Due to technical reasons the colour of this Servalyt® 6-7 appears darker than in the past. This is without any impact on quality and performance.		

⁽¹⁾ Gelelektrophorese wird auf Standard-PAG-Gelen durchgeführt.
 Gel electrophoresis is performed using a standard-PAG-gel.
 (5% T, 3% C, 4% Servalyt®, 0.1% (w/v) APS, 0.078% (w/v) TEMED,
 standard electrode solution # 42984 and # 42986)

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 26.05.2020

Christian Monsler

Marina Konter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.