

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Applikatorstreifen-Kit	Kat.Nr./Cat.No. :	42937
Applicator strips kit	Lot/Contr.No. :	161149

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Inhalt content		enthält je 1x Kat.Nr. 42989, 42914, 42915, 42899 contains 1 each of cat.nos. 42898, 42914, 42915, 42899	
Anwendung application		für Isoelektrische Fokussierung for Isoelectric focussing	
Mindesthaltbarkeit minimum shelf life			24.09.2018
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Datum/date: 13.12.2016

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.