

**A N A L Y S E N Z E R T I F I K A T**  
 CERTIFICATE OF ANALYSIS

<b>Applikatorstreifen-Kit</b> Applicator strips kit	<b>Kat.Nr./Cat.No. :</b> 42937 <b>Lot/Contr.No. :</b> 190884
--	---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Inhalt</b> content		enthält je 1x Kat.Nr. 42989, 42914, 42915, 42899 contains 1 each of cat.nos. 42989, 42914, 42915, 42899	
<b>Anwendung</b> application		für Isoelektrische Fokussierung for Isoelectric focussing	
<b>Mindesthaltbarkeit</b> minimum shelf life			18.01.2022
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
 We do not guarantee that the product can be used for a special application.  
 This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 05.09.2019

Christian Monsler

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.