

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Applikatorstreifen-Kit Applicator strips kit	Kat.Nr./Cat.No. : 42937 Lot/Contr.No. : 201174
--	---

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Inhalt content		enthält je 1x Kat.Nr. 42989, 42914, 42915, 42899 contains 1 each of cat.nos. 42989, 42914, 42915, 42899	
Anwendung application		für Isoelektrische Fokussierung for Isoelectric focussing	
Mindesthaltbarkeit minimum shelf life			21.08.2023
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 03.11.2020

Christian Monsler

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.