

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

| | | |
|-------------------------------|--------------------------|---------------|
| Applikatorstreifen-Kit | Kat.Nr./Cat.No. : | 42937 |
| Applicator strips kit | Lot/Contr.No. : | 211001 |

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|---|--------------------------|--|---------------------------|
| Inhalt content | | enthält je 1x Kat.Nr. 42989, 42914, 42915, 42899 contains 1 each of cat.nos. 42989, 42914, 42915, 42899 | |
| Anwendung application | | für Isoelektrische Fokussierung for Isoelectric focussing | |
| Mindesthaltbarkeit minimum shelf life | | | 02.04.2024 |
| Lagerung (°C) storage | | | +15 bis/ to +30 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 25.08.2021

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.