

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Elektrodendochte lang	Kat.Nr./Cat.No. :	42987
Electrode wicks long size	Lot/Contr.No. :	120820

Parameter parameter	Spezifikation specification	Ergebnis result
Material material	Filterkarton filter cardboard	entspricht corresponds
Länge length	240 x 6 x 1 mm	entspricht corresponds
Anwendung application	für Isoelektrische Fokussierung for Isoelectric focussing	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		07/ 2017
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 30.06.2014

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.