

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Elektrodendochte lang	Kat.Nr./Cat.No. :	42987
Electrode wicks long size	Lot/Contr.No. :	150643

Parameter parameter	Spezifikation specification	Ergebnis result
Material material	Filterkarton filter cardboard	entspricht corresponds
Länge length	240 x 6 x 1 mm	entspricht corresponds
Anwendung application	für Isoelektrische Fokussierung for Isoelectric focussing	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		01.06.2020
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 20.11.2019

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.