

**A N A L Y S E N Z E R T I F I K A T**  
CERTIFICATE OF ANALYSIS

<b>GEL-FIX® for PAG</b>	<b>Kat.Nr./Cat.No. : 42993</b>
<b>265 x 125 mm</b>	<b>Lot/Contr.No. : P140561</b>

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Beschreibung</b> description	klare, farblose Haftfolien aus Polyester, auf beiden Seiten aktiviert mit einer Acrylatschicht clear, colorless polyester sheets, activated on both sides to bind polyacrylamide gels		
<b>Größe (mm)</b> size	265 x 125		entspricht corresponds
<b>Dicke (mm)</b> thickness	0.18		entspricht corresponds
<b>Anwendung</b> application	zum Herstellen von Polyacrylamidgelen for producing polyacrylamide gels		
<b>Mindesthaltbarkeit</b> minimum shelf life	26.11.2017		
<b>Lagerung (°C)</b> storage	+15 bis/ to +30		

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b>  <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b>          We do not guarantee that the product can be used for a special application.          This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 17.01.2017**

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.