

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

SERVAGel(TM) IEF Starter Kit

Kat.Nr./Cat.No. : 43205

Lot/Contr.No. : P160255

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	Für die isoelektrische Fokussierung von Proteinen For the isoelectric focusing of proteins	entspricht corresponds
Inhalt content	SERVAGel(TM) IEF 3 – 10 SERVA IEF Anode Buffer, powder SERVA IEF Cathode Buffer, powder IEF sample buffer (2x) IEF Marker 3-10, Liquid Mix SERVA Violet 17	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		14.01.2017
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 12.07.2016

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.