

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

SERVAGel(TM) IEF Starter Kit	Kat.Nr./Cat.No. : 43205
	Lot/Contr.No. : P170105

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	Für die isoelektrische Fokussierung von Proteinen For the isoelectric focusing of proteins	entspricht corresponds
Inhalt content	SERVAGel(TM) IEF 3 – 10 SERVA IEF Anode Buffer, powder SERVA IEF Cathode Buffer, powder IEF sample buffer (2x) IEF Marker 3-10, Liquid Mix SERVA Violet 17	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		07.07.2017
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 08.03.2017

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.