

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>SERVAGel(TM) IEF Starter Kit</b>	<b>Kat.Nr./Cat.No. :</b>	<b>43205</b>
	<b>Lot/Contr.No. :</b>	<b>P250230</b>

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Anwendung</b> application	Für die isoelektrische Fokussierung von Proteinen For the isoelectric focusing of proteins	entspricht corresponds
<b>Inhalt</b> content	SERVAGel(TM) IEF 3 – 10 SERVA IEF Anode Buffer, powder SERVA IEF Cathode Buffer, powder IEF sample buffer (2x) IEF Marker 3-10, Liquid Mix SERVA Violet 17	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life		10.03.2026
<b>Lagerung (°C)</b> storage		+2 bis/ to +8

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b> <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b> We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
---

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 22.07.2025**

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.