

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

2D HPE(TM) Large Gel NF 12.5% Kit	Kat.Nr./Cat.No. :	43304
	Lot/Contr.No. :	141114 - 140860

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	Kit enthält 4 an Kunststoffolie gebunde Gele, einschließlich Lauf- und Äquibrierungspuffer, FS Wicks und Kühlkontaktflüssigkeit. Format: 255 x 200 mm, Dicke: 0.65 mm		
	Kit comprising 4 plastic-backed gels, including running and equilibration buffers, FS wicks and cooling contact fluid. Size: 255 x 200 x 0.65 mm		
Anwendungstest application test	SDS-Elektrophorese SDS-Electrophoresis	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			06/ 2016
Lagerung (°C) storage	Buffer Kit for 2D HPE(TM) Gels		+15 bis/ to +30
Lagerung (°C) storage	2D HPE(TM) Large Gel NF-12.5 %		+2 bis/ to +8

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 24.11.2014

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.