

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

---

<b>2D HPE(TM) Large Gel NF 12.5% Kit</b>	<b>Kat.Nr./Cat.No. :</b>	<b>43304</b>
	<b>Lot/Contr.No. :</b>	<b>141169-141180</b>

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Beschreibung</b> description	Kit enthält 4 an Kunststoffolie gebunde Gele, einschließlich Lauf- und Äquilibrationspuffer, FS Wicks und Kühlkontaktflüssigkeit. Format: 255 x 200 mm, Dicke: 0.65 mm		
	Kit comprising 4 plastic-backed gels, including running and equilibration buffers, FS wicks and cooling contact fluid. Size: 255 x 200 x 0.65 mm		
<b>Anwendungstest</b> application test	SDS-Elektrophorese SDS-Electrophoresis	entspricht corresponds	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			07/ 2016
<b>Lagerung (°C)</b> storage	Buffer Kit for 2D HPE(TM) Gels		+15 bis/ to +30
<b>Lagerung (°C)</b> storage	2D HPE(TM) Large Gel NF-12.5 %		+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

---

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 15.01.2015**

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.